

DAY 1

# Foundations & Amplicon Metagenomics

What is metagenomics • Types of sequencing • 16S rRNA gene  
Databases • Amplicon analysis pipeline

# Today's Learning Objectives

---

## 01 Define metagenomics

Explain the concept of studying genetic material from environmental samples without culturing.

## 02 Compare sequencing strategies

Distinguish amplicon from shotgun metagenomics and choose the right approach.

## 03 Understand the 16S rRNA gene

Describe its structure, variable regions, and why it is used as a phylogenetic marker.

## 04 Navigate key databases

Identify SILVA, Greengenes2, NCBI RefSeq and their roles in taxonomy assignment.

## 05 Trace the amplicon pipeline

Follow raw reads from QC to taxonomy using QIIME2, DADA2 and EMU.

## 06 Hands-on with real data

Apply FASTQC, SEQKIT and QIIME2 to human gut and soil 16S datasets.

# What is Metagenomics?

*"The study of genetic material recovered directly from environmental samples, without prior cultivation of organisms."*



## Culture-independent

Study microbes that cannot be grown in the lab — >99% of environmental microbes are unculturable.



## Community-level

Analyse entire microbial communities simultaneously: bacteria, archaea, viruses, fungi in one experiment.



## High-throughput

Next-generation sequencing generates millions of reads, enabling deep coverage of complex communities.

# Why Metagenomics? Applications Across Fields

## Human Health

Gut microbiome & disease (IBD, obesity, T2D)  
Probiotics & precision medicine  
Hospital-acquired infection surveillance

## Environmental Science

Soil health & agriculture  
Ocean & freshwater ecology  
Biodiversity monitoring & conservation

## Biotechnology

Novel enzyme discovery (bioprospecting)  
Bioactive compound mining  
Bioremediation of pollutants

## Food & Fermentation

Fermented food safety (dairy, beer, sauerkraut)  
Food spoilage detection  
Probiotic strain characterisation

## Agriculture

Plant-growth-promoting bacteria  
Soil carbon cycling & nutrient availability  
Crop protection & disease suppression

## Epidemiology

Outbreak surveillance & pathogen detection  
Antibiotic resistance gene tracking  
Zoonotic spillover monitoring

# Approaches in Metagenomics

## Amplicon Metagenomics

(Targeted)

PCR-amplify a marker gene (16S, 18S, ITS) from the community, then sequence the amplicons.

### ✓ Advantages

- Low cost per sample
- Deep coverage of marker gene
- Well-established pipelines
- Handles complex communities

### ✗ Limitations

- Only reveals who is there (taxonomy)
- PCR bias and chimera formation
- No functional information
- Primer mismatches miss taxa

THIS COURSE FOCUS →

## Shotgun Metagenomics

(Whole-genome)

Randomly fragment and sequence all DNA in the sample, no amplification step.

### ✓ Advantages

- Taxonomic + functional data
- No PCR bias
- Detect novel organisms
- Reconstruct genomes (MAGs)

### ✗ Limitations

- High cost & computation
- Host DNA contamination
- Requires large databases
- Rare taxa easily missed

## Metatranscriptomics

(RNA-based)

Sequence RNA (mRNA) from the community to capture active gene expression.

### ✓ Advantages

- Shows what is actively expressed
- Captures community response
- Reveals metabolic activity

### ✗ Limitations

- RNA unstable — fast processing needed
- Complex analysis
- rRNA dominates (~95%)

SECTION 2

# The 16S rRNA Gene

Structure • Variable Regions • Why It Matters

# The 16S rRNA Gene

The 16S rRNA gene encodes the small subunit of the prokaryotic ribosome. It is ~1,540 bp long and contains:

## 9 Conserved Regions

Shared across all prokaryotes — used for universal primer design

## 9 Variable Regions (V1–V9)

Evolve at different rates — provide species/genus-level resolution

## ~1,540 bp total

Full-length sequenced with PacBio/ONT; V3-V4 (~460 bp) most common with Illumina

## Universal marker

Present in all bacteria and archaea; absent from eukaryotes (use 18S/ITS instead)

## 16S rRNA Gene Variable Regions



| Region      | Forward Primer                | Reverse Primer                    | Read Length       |
|-------------|-------------------------------|-----------------------------------|-------------------|
| V4          | 515F<br>(GTGYCAGCMGCCGCGGTAA) | 806R<br>(GGACTACNVGGGTWTCTAAT)    | Illumina<br>2×250 |
| V3-V4       | 341F<br>(CCTACGGGNGGCWGCAG)   | 806R<br>(GGACTACNVGGGTWTCTAAT)    | Illumina<br>2×300 |
| V1-V3       | 27F<br>(AGAGTTTGATCMTGGCTCAG) | 518R (ATTACGCGGCTGCTGG)           | Illumina<br>2×300 |
| Full-length | 27F<br>(AGAGTTTGATCMTGGCTCAG) | 1492R<br>(TACGGYTACCTTGTTACGACTT) | ONT / PacBio      |

*Primer choice determines which microbial groups are detected and at what taxonomic resolution.*

# Sequencing Platforms for Metagenomics

## Illumina (Short Reads)

### Specifications

- Read length: 2×150 or 2×300 bp
- Throughput: up to 3 Tb/run (NovaSeq)
- Error rate: <0.1% (very accurate)
- Cost: ~\$50–200 per sample

### Common Uses

- V4/V3-V4 amplicon sequencing
- Shotgun WGS (most common)
- Large cohort studies

Tools: FASTQC, Trimmomatic, fastp, QIIME2, DADA2

## Oxford Nanopore (ONT)

### Specifications

- Read length: up to >1 Mb (avg 5–50 kb)
- Throughput: 50–300 Gb/run (PromethION)
- Error rate: ~5% (improving with R10 chemistry)
- Cost: ~\$500–800/flow cell (~\$50–150/sample multiplexed)

### Common Uses

- Full-length 16S sequencing
- Rapid pathogen detection
- Direct RNA sequencing

Tools: wf-16S (EPI2ME), EMU, guppy/dorado, minimap2

## PacBio HiFi

### Specifications

- Read length: 10–25 kb (HiFi)
- Accuracy: >99.9% (CCS mode)
- Lower throughput than Illumina
- Higher cost per Gb

### Common Uses

- Full-length 16S/18S (near-perfect accuracy)
- Genome assembly (hybrid approach)
- Large structural variation

Tools: pbmm2, DADA2 (with HiFi reads)



# Amplicon vs Shotgun vs Long-Read: When to Use Each

|                   | Amplicon (Short-Read)              | Shotgun WGS (Short-Read)            | Long-Read (ONT / PacBio)                   |
|-------------------|------------------------------------|-------------------------------------|--------------------------------------------|
| What is sequenced | 16S / ITS / 18S amplicon           | All DNA in sample (WGS)             | Full-length amplicon or WGS                |
| Cost per sample   | ~\$50–100 (cheap)                  | ~\$100–300 (moderate–high)          | ~\$50–150 ONT (falling fast)               |
| Resolution        | Genus/species (V4); species (FL)   | Strain-level; functional genes      | Species-level (FL 16S); strain (WGS)       |
| Error profile     | <0.1% substitution (Illumina)      | <0.1% substitution (Illumina)       | ~5% indel-rich (ONT); <0.01% (PacBio HiFi) |
| Functional info   | None (marker gene only)            | Full metabolic potential + AMR      | Full (WGS); none (amplicon)                |
| Best for          | Community profiling, large cohorts | Assembly, MAGs, functional analysis | Species ID in clinical; complex assembly   |
| Key tools         | QIIME2, DADA2, EMU                 | Kraken2, MEGAHIT, MetaBAT2          | Dorado, EMU, wf-16S, metaFlye              |

# Primer Choice for Amplicon Sequencing

## 16S rRNA — V4

V4 (~253 bp) · 515F / 806R

✓ Most widely used ✓ Illumina 2×150/250 ✓ EMP standard ✓ Good archaea coverage

✗ Misses some Thaumarchaeota ✗ Low resolution below genus level

## 16S rRNA — V1-V3

V1-V3 (~500 bp) · 27F / 518R

✓ This course's dataset (PRJNA683885) ✓ Longer amplicon → more resolution ✓ Classic legacy marker

✗ No standard pre-trained SILVA classifier — must build locally ✗ Requires 2×300 bp reads

## 16S rRNA — V3-V4

V3-V4 (~460 bp) · 341F / 806R

✓ Better taxonomic resolution ✓ Common in human microbiome studies ✓ SILVA v138.2 classifiers

✗ Requires 2×300 bp reads ✗ More PCR cycles → more chimeras

## 16S rRNA — Full-length

V1-V9 (~1500 bp)

27F / 1492R

✓ Species-level resolution ✓ ONT & PacBio HiFi only ✓ EMU & wf-16S designed for this

✗ Long reads required (ONT/PacBio) ✗ Higher error rate (ONT) ✗ Cannot use Illumina 2×300

## ITS (Fungi)

ITS1 or ITS2 (~300 bp)

ITS1F / ITS2 or ITS3F / ITS4

✓ Standard for mycobiome studies ✓ UNITE v10.0 database ✓ High species-level resolution

✗ ITS length variability complicates DADA2 ✗ Consider ITSxpress for trimming

# Sequencing Error Profiles & Quality Control

## Illumina (MiSeq/NovaSeq)

Error type: Substitution errors

Quality drops at 3' end; systematic run-specific errors

QC/handling: DADA2 error model per run; fastp Q20 filter

<0.1%

## Oxford Nanopore (R10)

Error type: Indel errors (ins/del)

Homopolymer runs (poly-A/T stretches) are hardest; R10 much better than R9

QC/handling: EMU / minimap2 handle indels; Dorado recalibration

~5%

## PacBio HiFi (CCS)

Error type: Substitution errors

Multiple passes produce near-perfect consensus; lowest error of any long-read

QC/handling: DADA2 with HiFi-specific parameters; pbmm2

<0.01%

SECTION 3

# Reference Databases

SILVA • Greengenes2 • NCBI RefSeq • UNITE • GTDB

# Key Reference Databases

## SILVA

v138.2

### 16S / 18S / 23S rRNA

Gold standard rRNA database (July 2024).  
Curated alignments for all three domains of life.  
Over 9.4 million SSU sequences. QIIME2 classifiers trained on SILVA are widely used.

## Greengenes2

2024.09

### 16S rRNA

Rebuilt integrating Wol2 (Web of Life 2), a phylogenomic reference tree from ~10,000+ genomes.  
Better species-level resolution. Default for recent QIIME2 tutorials.

## NCBI RefSeq

Live

### Full genomes

Comprehensive repository of reference sequences.  
Used for Kraken2 database, BLAST searches, and downloading reference genomes.

## UNITE

10.0

### ITS (Fungi)

Species Hypothesis (SH) based fungal database. Essential for amplicon studies targeting ITS1/ITS2 regions for mycobiome analysis.

## GTDB

r220

### Whole genomes

Genome Taxonomy Database, phylogenomics-based prokaryote taxonomy (April 2024 release). Standardised nomenclature used with GTDB-Tk for MAG classification.

## eggNOG 5.0

5.0.2

### Proteins / COG

Hierarchical orthologous groups. Used by eggNOG-mapper for rapid functional annotation. Maps genes to COG categories, GO terms and KEGG pathways.

SECTION 4

# The Amplicon Analysis Pipeline

From raw reads to taxonomy and feature tables

# Amplicon Analysis Workflow

## 01 Raw Reads

### FASTQC + MultiQC

Quality control of raw FASTQ files — check quality scores, adapter content, GC bias.

## 02 Trimming & Filtering

### Trimmomatic / fastp

Remove low-quality bases, adapters; filter very short reads.

## 03 Import to QIIME2

### qiime tools import

Convert FASTQ files to QIIME2 artifact format (.qza). Define manifest file.

## 04 Denoise & Chimera Removal

### DADA2 / Deblur

Error correction and chimera removal; generate ASVs (Amplicon Sequence Variants). DADA2 models sequencing errors per run.

## 05 Taxonomy Assignment

### q2-feature-classifier

Classify ASVs using SILVA or Greengenes2 reference database.

## 06 Diversity & Visualisation

### q2-diversity / R

Alpha/beta diversity, taxonomy bar plots, statistical testing.



### Long-Read Track (ONT):

Basecall (Dorado) → Full-length 16S reads → EMU (species-level abundance) or wf-16S (EPI2ME Nextflow workflow) | Skips QIIME2 DADA2 — uses alignment-based error handling instead

# Oxford Nanopore 16S Tools

For long-read ONT 16S sequencing, different pipelines are used compared to Illumina:

## EMU

*Exact Mapping Utilities*

Species-level 16S abundance estimation for ONT reads. Uses multi-mapping alignment to handle errors. Supports full-length 16S.

### Pipeline steps:

1. Input: raw FASTQ from ONT
2. Maps reads to EMU default database (SILVA-based)
3. Outputs: species-level abundance table
4. Compatible with QIIME2 via q2-emu plugin

## wf-16S

*EPI2ME Nextflow workflow*

Official Oxford Nanopore Technologies workflow for 16S/18S amplicon analysis. Nextflow-based, portable, cloud-ready.

### Pipeline steps:

1. Input: FASTQ or POD5 files
2. Basecalling → Demultiplexing → Classification
3. Uses Minimap2 + SILVA/Greengenes2
4. Outputs: abundance table + report

## Dorado

*Basecaller (ONT)*

GPU-accelerated basecaller from ONT. Converts raw signal (POD5) to FASTQ. Supports demultiplexing and modified base detection.

### Pipeline steps:

1. Input: POD5 raw signal files
2. Basecalling model selection (fast/hac/sup)
3. Output: FASTQ with quality scores
4. Replaces older Guppy basecaller



# Tools We Use in This Course

---

## Day 1–2

Amplicon & Diversity

FASTQC · MultiQC · SEQKIT · QIIME2 · DADA2 · EMU · epi2me/wf-16S · phyloseq · vegan · ggplot2

## Day 3

Shotgun & Assembly

fastp · SEQKIT · Bowtie2 · BWA · Samtools · Kraken2 · Bracken · MetaPhlAn4 · MEGAHIT · SPAdes · QUAST

## Day 4

MAGs

MetaBAT2 · DAS\_Tool · CheckM2 · GTDB-Tk · Bowtie2 · Samtools

## Day 5

Functional Analysis

Prodigal · Prokka · eggNOG-mapper · HUMAnN3 · R (pathway plots)

*All tools available via module load. See `environment/metagenomics.yml` and `environment/qiime2.yml` for installation.*

# Today's Hands-On Session

## 1 QC Your Reads

Run FASTQC on raw gut microbiome 16S reads. Generate a MultiQC summary report. Inspect quality plots and identify potential issues.

## 2 Explore with SEQKIT

Use seqkit stats to get read count, length distributions, and GC content. Practice subsetting and filtering with seqkit grep, head, sample.

## 3 QIIME2 Amplicon Workflow

Import paired-end reads → run DADA2 denoising → assign taxonomy with SILVA classifier → generate taxonomy bar plots and feature tables.

## 4 ONT 16S with EMU (Bonus)

If time allows, process a small ONT 16S dataset using EMU. Compare taxonomy profiles to the Illumina-based results from step 3.

# Day 1 Summary

---

## Metagenomics

Culture-independent study of genetic material from complex environmental communities.

## Amplicon sequencing

PCR-amplifies a marker gene (16S rRNA) for community profiling — cost-effective and well-established.

## 16S rRNA gene

~1,540 bp; 9 variable regions; V4 most common for Illumina; full-length for ONT/PacBio.

## Databases

SILVA & Greengenes2 for 16S taxonomy; NCBI for reference genomes; GTDB for standardised taxonomy.

## Pipeline

FASTQC → QIIME2 import → DADA2 (ASVs) → taxonomy classification → diversity analysis.

## ONT 16S

Full-length 16S with EMU or wf-16S provides species-level resolution not achievable with short reads.

# KEY REFERENCES

## **Bolyen et al. (2019)**

Reproducible, interactive, scalable microbiome data science using QIIME 2

*Nature Biotechnology* 37: 852–857

## **Callahan et al. (2016)**

DADA2: High-resolution sample inference from Illumina amplicon data

*Nature Methods* 13: 581–583

## **Quast et al. (2013)**

The SILVA ribosomal RNA gene database project

*Nucleic Acids Research* 41: D590–D596

## **McDonald et al. (2024)**

Greengenes2 unifies microbial data in a single reference tree

*Nature Biotechnology* 41: 1730–1737

## **Curry et al. (2022)**

Emu: species-level profiling of full-length 16S rRNA ONT data

*Nature Methods* 19: 845–853

## **Shen et al. (2016)**

SeqKit: a cross-platform ultrafast toolkit for FASTA/Q manipulation

*PLOS ONE* 11: e0163962